

THERAPEUTICS DISTRIBUTORS REBEL FOOT FLEX INSTRUCTION GUIDE

Step 1: Tension Strap is placed into Sock Channel Guide.



Step 2: Press the Channel Guide Velcro closed.



Step 3: Put Sock on foot.



Step 4: With knee hole to centered on kneecap, wrap Lower Adjustment Strap around front of leg; pull and tighten.



Step 5: Fasten Lower Adjustment Strap.



Step 6: Wrap the Upper Adjustment Strap around the back of the leg to the front; pull, tighten, then fasten.



Step 7: Attach Tension Strap buckle to Knee Assembly (receptor) buckle on left side.



Step 8: Attach Tension Strap buckle to Knee Assembly (receptor) buckle on right side. *Note: Letters should be oriented as shown.*



Step 9: Place Tension Strap under the Alignment Tab.



Step 10: Open the Cam Lock then pull up on Tension strap.



Step 11: Close the Cam Lock to secure desired tension; use Velcro on Tension Strap end to secure excess strap to Knee Assembly.



Step 12: The Foot Flex demonstrated when instructions are completed.

