

THERAPEUTICS DISTRIBUTORS REBEL FOOT FLEX INSTRUCTION GUIDE

Step 1: Place Tension Strap into Sock or Flex Slipper Channel Guide. *Note: Keep Tension Strap letters facing outward.*



Step 2: Press the Channel Guide Velcro closed on bottom of Sock or Flex Slipper.



Step 3: Put Sock or Flex Slipper on foot. *Note: If your Foot Flex is equipped with a Flex Slipper, see the additional instruction page.*



Step 4: Center Knee hole on kneecap; Wrap Lower Adjustment Strap around front of leg; Pull to tighten.




Step 5: Fasten Lower Adjustment Strap.



Step 6: Wrap the Upper Adjustment Strap around the back of the leg to the front; pull to tighten, then fasten.



Step 7: Attach Tension Strap buckle to Knee Assembly (receptor) buckle on LEFT side.



Step 8: Open the Cam Lock. *Note: Keep letters on Tension Strap oriented to right side so they are visible.*



Step 9: Attach Tension Strap Buckle to Knee Assembly (receptor) Buckle on RIGHT side.



Step 10: Place Tension Strap Under the Alignment Tab.



Step 11: Close the Cam Lock to secure desired tension; use Velcro on Tension Strap end to secure excess strap to Knee Assembly.



Step 12: The Foot Flex demonstrated when instructions are completed.

